

-30/+10°C

REMEMBER TO GET A QUOTE ON OPTIONS

AVAILABLE IN STAINLESS STEEL EXTERIOR



LABORATORY REFRIGERATORS AND FREEZERS

# LFF series (Combi Range)

This model combines all the technology and features of the LR/LF Series, but it is a real space saver, since it is equipped with both refrigerator and freezer, in one single unit. Each compartment has its own controller of course, in order to program them according to your own requirements and specifications.

All the models of the LR/LF range are re-enforced with extra holes on the top panels to enable the unit to cool down even faster.



LFF 270



### LFF series

- White or stainless steel exterior
- Stainless steel interior
- Automatic door closing
- Air forced cooling
- Adjustable shelves
- Castors
- Door key lock
- Reversible door

### Controller features

- Micro processor controller with digital display
- Approx. 72 hours battery back up for alarms, loggings and temperature display in case of power cut
- Visual and acoustic alarm
- Adjustable high/low temperature alarm
- Power failure alarm
- Probe failure alarm
- Instrument failure alarm
- Open door alarm
- Contact for remote alarm
- Prepared for GSM alarm
- Integrated data logger (software included)
- RS485/232 Interface
- Computer USB data read out
- Auto cycle if probe failure
- Ambient temperature display
- Strows all alarms
- Direct download/upload on/ from USB memory stick
- 3-level password protected
- Battery level indication
- Integrated memory

Scan me and see more about Combi Range



www.arctiko.com

Model	LFF 270**	LFF 660**
Item code (DAI no.)	DAI 0235/0236*	DAI 0265/0266*
Dimensions exterior WxDxH (mm)	520x700x1997	720x860x1997
Net capacity (L)	160 / 160	290 / 290
Temperature range (°C)	+1/+10 / -10/-30	+1/+10 / -10/-30
Power supply (V)	230	230
Frequency (Hz)	50	50

We reserve the right to change specifications without notice. Subject to confirmation and availability. \* DAI number for unit with stainless steel exterior. \*\* For stainless steel exterior please add -ST to the existing product name.